


PURPOSE T Pressure Ulcer Risk Assessment

Mattress Selection Guide at DBTH

Skin Integrity Team



The Mercury Advance Hybrid Mattress provides a solution to 'step up' to that of an air mattress when clinically required and 'stepped down' to a foam mattress as the patient's condition improves, whilst using the same mattress. Here is a user video:
<https://www.youtube.com/watch?v=h1a19Ufoxds>

The Quattro Plus Dynamic mattress is a high specification air mattress to aid with the management of up to a category 4 pressure ulcer. Here are some user video for the Quattro Plus Dynamic:
<https://www.youtube.com/watch?v=h1a19Ufoxds>
<https://www.youtube.com/watch?v=tOb4NJcWfFI>
<https://www.youtube.com/watch?v=BHR1ZJ2OHSU>